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PROSPDEFEXP0123-rev011

IMPORTANT INFORMATION. CAREFULLY READ THIS DOCUMENT BEFORE USING THE PRODUCT.

INSTRUCTIONS FOR USE

Phibo® instruments and prosthetic components.

Pilar ProUnic Plus® Phibo® Avantblast® TSA® Advance TSA® BNT® TSH® Phibo® Duplit® Softtissue® International Phibo Group® IFIC® VK® Bonetech® General® Esthetic Tissue® Phibo Esthetics® Phibo Surgical® Phibo Prostodontics® Phibo Scientific® Dual-Press® Non Stop® TSP® Aurea® Aurea® Evo are trademarks of Phibo Dental Solutions, S.L.

This document contains information for the use of prosthetic components and instruments of the Phibo® dental implant systems: TSA®, TSA® Advance, TSH®, BNT®, Aurea®, Aurea®, Evo.

The components of Phibo® dental implant systems should be used only by appropriately trained professionals. For detailed information about product specifications for use please refer to the respective surgical and prosthetic procedures.

INFORMATION ON LIABILITY AND WARRANTY

The user should ensure that the product employed is suitable for its intended use, particularly in case of procedures not explicitly recommended. The user may not disclaim this liability based on those instructions for use or procedures. The products in Phibo® dental implant systems should only be used with original components and instruments in accordance with the corresponding instructions and recommendations. The use of non-original Phibo® products, instruments or components that come into contact with the ones mentioned in the Phibo® catalog and procedures shall automatically render warranties on products manufactured by Phibo Dental Solutions, S.L. void.

The professional in charge of clinical treatment should guarantee that original components be used throughout the clinical and laboratory process. To obtain any kind of warranty, the professional in charge of the clinical treatment and of the patient should provide all the information required by Phibo® to that end.

1. DESCRIPTION

PROSTHODONTIC COMPONENTS AND INSTRUMENTS

The Phibo® equipment is made of stainless steel, tungsten carbide and titanium. Surgical drills are designed to be used a total of ten times. Inadequate maintenance or lack of cleaning and disinfection may reduce the number of uses, apart from causing treatment failure. All instruments are color-coded or laser-marked for easy identification and appropriate use. Using instruments of other commercial brands implies the loss of any kind of guarantee on the implants and the rest of the products. Phibo® prosthetic components are manufactured of titanium and biocompatible plastic material. Some Phibo® abutments and screws are color coded for their easy identification and appropriate use according to each implant series.

2. INDICATIONS

INSTRUMENTS AND ATTACHMENTS

Phibo® instruments are designed to support care planning. Osseous bed preparation and Phibo® dental implants insertion in the patient's maxilla or mandible. The instruments include some screwdrivers and instruments for prostheses and orthodontics components, as per the product catalog. Phibo® abutments are fixed to Phibo® dental implants to provide support to the prosthesis. Prostheses may be single, partial or total, and may be screwed on or cemented to the abutment, or may be removable.

Plastic prosthetic components for provisional use are designed as a support for temporary restorations for periods of no longer than 60 days. These products should be used with the appropriate Phibo® dental implants and abutments in accordance with the corresponding surgical and prosthetic procedures.

Abutments with an embrioging handle can be removed by hand.

It is important to regulate both the counter-angle in case of mechanical insertion, and the torque wrench in case of manual insertion, proceed as indicated in the surgical procedure.

Inserting the implant and exceeding that predetermined force will be a sufficient indication of the need to carry out the complete, and not just partial, surgical sequence again as defined in the surgical procedure description.

PROSTHODONTIC REHABILITATION

Prosthetic procedures describe in detail the precautions to be taken during therapy. The prosthesis and rehabilitation designs should always be based on the type of prosthetic rehabilitation required by the patient and which should ensure therapeutic success and patient expectations.

6. ADVERSE EFFECTS

Implantology include side effects that are documented in the specialized scientific literature published in the dentistry sector. However, the most relevant effects are: post-operative discomfort, local inflammation, local or systemic infections, difficulty speaking, bone loss and fractures, implant loss, damage to adjacent teeth, fracture of implants and prosthetic components, and damage to dental nerves, among others.

3. CONTRAINDICATIONS

It is necessary to carry out a medical pre-operative examination of the patient in order to determine risk factors during the intervention, implant insertion or removal. Dental implants should not be used in patients who lack the medical clearance necessary to undergo implantological therapy and rehabilitation. In case of patients with localized or systemic risk factors that may impair the bone or soft tissue healing process, the professional in charge should assess the soft and potential risks of treatment.

RELATIVE: age, stress, smoking, pregnancy, bone deficiency, alcoholism, drug abuse and addictions in general, lack of oral hygiene, periodontal disease, among others.

ABSOLUTE: endocrine (decompensated diabetes mellitus, hyperthyroidism), blood dyscrasias contraindicating surgical treatments, cardiovascular and/or end-stage pathologies, infectious diseases, radiotherapy, corticotherapy and anticoagulant agents, epilepsy and psychological factors.

4. STORAGE AND HANDLING

Phibo products should be stored at a temperature between 10 and 40°C in a dry, clean place protected from adverse conditions.

5. WARNING

Treatment planning and dental implant positioning require specific odontological training. It is recommended that clinical users attend hands-on training courses to learn appropriate techniques, including treatment-associated biomechanical and radiographic requirements. Before positioning Phibo® dental implants or their prosthetic elements it is necessary to be familiar with the corresponding surgical and prosthetic procedures.

The patient should have the bone volume and bone quality adequate for insertion of the necessary implants and for support of the functional loads foreseen in the service. The professional in charge of the implantological therapy should, by means of correct rehabilitation planning, ensure an adequate safety margin, together with teeth and vital structures. Otherwise, serious damage may be caused to vital anatomical structures, with temporary and/or permanent lesions, as well as harm to the patient's health.

Each dental implant system has its own design characteristics encompassing implants, prosthetic components and instruments. The use of inappropriate components or components made by other manufacturers may result in component mechanical failure, tissue damage or definite aesthetic results, due to lack of compatibility with specifications Microdesign & Macrodesign.

The procedure for employing drills, counter-angle bone taps and other instruments needed for implant placement are detailed in the respective surgical procedure descriptions.

Implant placement and prosthetic planning should adapt to the patient's individual conditions, particularly with respect to the correct distribution of forces. A passive adjustment should be attained in prosthetic rehabilitation, with occlusion adjustment in the opposite maxillary bone, and excessive lateral forces should be avoided. An insufficient number of implants, a wrong size selection or an inappropriate positioning for bearing and transmitting the expected loads may result in mechanical failure of the implant, the abutment or the abutment screws due to overload or fatigue and the substantial loss of surrounding bone. The lack of adequate residual bone quantity and quality, onset of infection or diseases in general and changes in patient habits are some potential causes for the failure of osseointegration and therapy. Lack of bone or soft tissue may result in an unfavorable implant insertion and a less-than-desired esthetic outcome. Inadequate prosthetic rehabilitation may cause rehabilitation failure.

Re-use of disposable products may damage their technical specifications, and also carries the risk of tissue infection, surgical failure, prosthetic failure and/or injury to the patient's health.

6. DISINFECTION AND STERILIZATION

Phibo® instruments and prosthetic components are not supplied sterile. Do not use products whose packaging is damaged or previously open. Prosthetic components and instruments to be used in the mouth should be cleaned, disinfected and sterilized. Reusable instruments and components must also be cleaned, disinfected and sterilized after use.

IMMEDIATELY AFTER IMPLANTATION: Disassemble instruments comprising several parts into their components according to the instructions for use (e.g. wrench). Damaged or dull instruments should be set aside (sharp instruments should be replaced after a maximum of 10 uses) and disinfected and cleaned separately.

MANUAL CLEANING AND PREPARATION OF UNITS

Brush and rinse under running water to remove excess dirt in all zones and crevices of instruments for between 20 and 30 seconds until all residue is removed.

DISINFECTION OF UNITS

Immerse the instruments in an adequate disinfectant bath, strictly following the manufacturer's instructions regarding dose/concentration, immersing time and temperature.

Instrument Instrumental PRO (universal disinfectant for instruments), a disinfectant with CE marking especially indicated for dental health material is recommended.

UNIT DISINFECTION PROCESS

Disinfection may be manual or automatic. The instruments should not be in contact with one another.

1. Prepare the disinfectant solution using tap water and

following the manufacturer's instructions for use to achieve the final concentration desired (e.g. 5% Instrumental Instrumental PRD).

2. Mix the solution by stirring gently and then pour it into the ultrasonic bath tray.

3. Immerse the clean and thoroughly rinsed instruments in the bath for 30 minutes at 60°C with a power setting of 50 W (40 KHz), making sure that they are fully submerged and in contact with the disinfectant.

4. Rinse the material with sterile water for 30 seconds and dry thoroughly under aseptic conditions (with gauze, paper, filtered air, etc.). Do not use hydrogen peroxide, oxidizing acids (nitric acid, sulphuric acid, oxalic acid, etc.) or any product with a high chlorine content.

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