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### IPS e.max® CAD by Phibo®

### INSTRUCTIONS FOR USE

IMPORTANT INFORMATION. PLEASE READ THIS DOCUMENT CAREFULLY BEFORE USING THE PRODUCT.

### INSTRUCTIONS FOR USE

Phibo® is a registered trademark of Phibo Dental Solutions, S.L.  
IPS e.max CAD is a registered trademark of Ivoclar Vivadent AG.

This document contains information for the use of Phibo® prosthetic components.

## PRODUCT DESCRIPTION

IPS e.max CAD is a ceramic block for the CAD-CAM technique used in the restoration of individual teeth.

The following restorations are made with this material:

- Veneers.
- Inlays.
- Onlays.
- Partial crowns.
- Anterior and posterior crowns.
- Primary telescopic crowns.

## MATERIAL

Lithium disilicate glass-ceramic for the CAD-CAM technique.  
The IPS e.max CAD blocks are mainly composed of SiO<sub>2</sub>.  
Additional components: Other components: Li<sub>2</sub>O, K<sub>2</sub>O, MgO, Al<sub>2</sub>O<sub>3</sub>, P<sub>2</sub>O<sub>5</sub> and other oxides.

The restorations made with IPS e.max CAD are designed to be adapted to each individual patient, and as such are single-use products. Reuse of single-use products may cause their deterioration, at the risk of tissue infection, surgical or prosthodontic failure and/or deterioration of patient health.

## CONTRAINDICATIONS

- Complete veneers on posterior crowns.
- Very deep subgingival preparations.
- Patients with very few remaining teeth.
- Bruxism.
- Any other use not specified in the indications.

## SIDE EFFECTS

The product should not be used to prepare restorations if the patient is allergic to any of the components of IPS e.max.

## RESPONSIBILITIES AND WARRANTY

These prosthetic components are part of a global concept and should only be used following the instructions and recommendations of Phibo®.

These prosthetic components should only be used by professionals with adequate training. In order for warranties to be valid, the clinician responsible for treatment must supply the pertinent information required by Phibo®.

Clinicians using Phibo® prosthetic components are responsible for establishing whether or not the product is indicated for each patient and each set of circumstances. Phibo® accepts no liability, either implicit or explicit, and will not be responsible for any damage, whether direct, indirect, punitive or of any other nature, produced by or arising in relation to errors of professional judgment or practice in the use of Phibo® products. Clinicians are also obliged to keep up to date with all the latest developments and applications regarding Phibo® products. In case of doubt, please contact Phibo®. Since the processing and surgical application of this product are under the control of the clinician, they are also the responsibility of the clinician. We accept no responsibility for damage arising from this cause.

## INSTRUCTIONS FOR USE FOR THE PROSTHETIC PROFESSIONAL

### 1. Production

An imprint-based model is prepared. Using the 3Shape system by Phibo®, the model is scanned and the IPS e.max by Phibo® restoration is designed. Once designed, it is sent to the Phibo® production center.

For optimum tone, the blocks are chosen based on the desired dental tone, the stump tone and/or the abutment tone.

### 2. Finish

If adjustments are needed, polish the IPS e.max CAD restorations while they are in the pre-crystallization state (blue color). Only adequate instruments should be used, operating at low rpm and low pressure in order to avoid chipping, and particularly fracture of the edges. Overheating of the glass-ceramic is to be avoided. Test the restorations on the stumps and complete finishing with great care. Check the proximal and occlusal points. Make sure that the minimum thicknesses are preserved even after making minor corrections. Design surface textures. Always clean the restoration in an ultrasound bath or with an air pistol before each crystallization. Do not sand the restorations with Al<sub>2</sub>O<sub>3</sub> or polishing pearls.

### 3. Crystallization and cosmetic/glassing procedure

Follow the instructions of Ivoclar Vivadent for the crystallization and cosmetic procedures according to the restoration involved and the processing technique used (cosmetic procedure, cut-back, stratification).

## WARNING

- Avoid inhaling the ceramic powder during the process: use aspiration equipment and a mask.
- The IPS Ceramic Etching Gel contains hydrofluoric acid. Always avoid contact with the skin, eyes and clothing, since this material is extremely toxic and corrosive. The gel has been developed for professional use only, and must not be used intraorally (in the oral cavity).
- IPS Contrast Spray Labside must not be used intraorally.

This leaflet was revised and approved in 2013-05-30. INSUSOEMAX\_rev003

Instructions for use  
IPS e.max® CAD  
by Phibo®

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IPS e.max CAD è un marchio registrato di Ivoclar Vivadent AG.