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
.....

### ENGLISH

### PHIBO CAD-CAM INTERMEDIATE STRUCTURES

IMPORTANT INFORMATION. READ THIS DOCUMENT CAREFULLY BEFORE USING THE DEVICE

### INSTRUCTIONS FOR USE

Phibo® is a registered trademark of Phibo Dental Solutions, S.L. This document contains information for the use of de Phibo® prosthetic structures.

### PRODUCT DESCRIPTION

Intermediate structure manufactured using CAD-DAM techniques for the creation of customized prostheses. These structures can be used for cemented and screwed restorations (single and multiple) and/or full contour (monolithic structures without porcelain veneering)-bars for overdentures for removable prostheses, and hybrid structures in different materials.

PMMA is a material used exclusively for provisional restorations and its period of use in the mouth is 6 months.

### MATERIALS

These structures can be manufactured in the following materials:

- Titanium.
- Gold titanium.
- Zirconia.
- Translucent Zirconia.
- Cobalt-Chrome.
- PMMA.

### PRECAUTIONS

Restorations produced using Phibo® CAD-CAM intermediate structures are custom-made for each patient; therefore, without any circumstances should they be used for another patient. Their re-use could imply a possible deterioration of their properties, which implies a risk of tissue infection, prosthodontic failure and/or deterioration of the patient's health.

### CONTRAINDICATIONS AND SIDE EFFECTS

Patients allergic to any of the materials described above. Any other use not specified in the product description.

### INFORMATION ON LIABILITY AND GUARANTEE

These prosthetic components are part of a global concept and should only be used in accordance with the instructions and recommendations of Phibo®. These prosthetic components should only be used by professionals with specialised training in dental prosthesis. In order for warranties to be valid, the clinician and technician responsible for the treatment shall provide all information required by Phibo®. The prosthetic laboratory is responsible for the design of the intermediate structure, as well as the finishing of the prostheses. Phibo® declines any liability, either implicit or explicit, and will not be held responsible for any direct, indirect, punitive or other type of damages caused by, or in relation to, any error in the judgement or practice of the professional in the use of Phibo® products. The clinician is responsible for determining whether a device is indicated or not for each patient and each situation. Phibo® is responsible for the manufacture of the intermediate structures designed by the prosthetic laboratory and prescribed by the clinician.

### PRECAUTIONS DURING THE PROCEDURE:

Close collaboration between the clinician and the dental technician is essential to ensure a successful comprehensive dental treatment. It is of special importance to achieve a correct distribution of forces through adaptation and passive adjustment of the prosthesis to ensure occlusion with the opposing jaw and preventing excessive transverse load, especially in cases of immediate load.

In the case of crowns or monolithic structures in Translucent Zirconia, special attention should be given to the design of the occlusal face.

### HANDLING PROCEDURE FOR THE DENTAL TECHNICIAN

#### 1. Finishing

- a) Procedure for finishing structures produced in Zirconia:
  - Check the adjustment of the piece on the model.
  - For products that require porcelain veneering, use materials that meet prosthetic criteria and follow the recommendations provided by the manufacturer of the ceramics, as well as the instructions given by Phibo® for the preparation of the structure and porcelain veneering on Zirconia.
  - For monolithic products in Translucent Zirconia, use suitable cosmetics to finish the prosthesis and follow the manufacturer's recommendations, as well as the Phibo® protocol for applying cosmetics to this type of structure.
  - In the event of rehabilitating screwed Zirconia crowns/ abutments with a Titanium base, make sure they are positioned correctly (there is only one correct position) and that the cement has polymerised completely before handling them. The fit of the prosthesis should be checked on the model before cementing the Titanium base to the crown or abutment.
- b) Procedure for finishing structures produced in Titanium:
  - Check the adjustment of the piece on the model.
  - Check the adjustment of the piece on the model.
  - c) Procedure for finishing structures produced in Cobalt-Chrome:
    - Check the adjustment of the piece on the model.
    - For products that require porcelain veneering, use materials that meet prosthetic criteria and follow the recommendations provided by the manufacturer of the ceramics, as well as the instructions given by Phibo® for the preparation of the structure and porcelain veneering on Cobalt-Chrome.
  - d) Procedure for finishing structures manufactured in PMMA:
    - Check the adjustment of the piece on the model.

#### 2. Cleanliness and disinfection.

- a) Manual cleaning and preparation of the structures: Brush all the surfaces and crevices of the device and rinse under running water for 20 to 30 seconds until all remains have been eliminated.
- b) Disinfecting structures

Immerse the device in a suitable disinfectant bath, strictly following the manufacturer's instructions: dosage/ concentration, time and temperature. The recommended disinfectant is Instrunet Instrumental PRD with CE mark (universal disinfectant for instruments), especially indicated for dental instruments and devices.

Disinfection process:
Disinfection may be manual or automatic. The different structures should not be in contact with each other.

  - Prepare the disinfectant solution using tap water and following the instructions provided by the manufacturer to obtain a final concentration of Instrunet Instrumental PRD of 5%.
  - Mix the solution by stirring gently and transfer it to the ultrasonic bath tray.
  - Immerse the perfectly cleaned and rinsed device in the bath for 30 minutes at 60°C and a power of 50 W (40 KHz), making sure it is completely submerged and in contact with the disinfectant.
  - Rinse the material with sterile water for 30 seconds and then dry carefully under aseptic conditions (gauze, paper, filtered air, etc.). Do not use hydrogen peroxide, oxidant acids (nitric acid, sulphuric acid, oxalic acid, etc.) or substances with high chlorine content.

#### 3. Structure placement

- In the event of **screwed implants**, make sure to screw the structure according to the table below showing the recommended torque for each type of material used:
  - Cobalt-Chrome on direct to implant connections: 35 Ncm
  - Cobalt-Chrome on connections on abutment: 25 Ncm
  - Titanium and Golden titanium on direct to implant connections: 35 Ncm
  - Titanium and Golden titanium on connections on abutment: 25 Ncm
  - Zirconia and Translucent Zirconia on direct to implant connections with Titanium base: 35 Ncm
  - Zirconia and Translucent Zirconia on direct to implant connections without Titanium base: 25 Ncm
  - Zirconia and Translucent Zirconia on connections on abutment: 25 Ncm
  - PMMA in all connections: 15 Ncm
- In the event of **screwed rehabilitations** on implants other than Phibo®, make sure to follow the recommendations of such manufacturer.
- It is essential to use the appropriate screw driver for each case. The available Phibo® screw drivers are: 1.0 mm hex driver, 1.25 mm hex driver, hexalobular (star) screw driver, Axis® screw driver for connections with correction of the screw angle, or Conhex screw driver.

### DESIGN AND PRODUCTION PARAMETERS

All Phibo® intermediate structures must comply with certain minimum design parameters. The table below indicates these parameters, depending on the type of material and rehabilitation to be performed.

Cemented space			
Name	Edge line adjustment (mm)	Angle adjustment (deg)	Wall thickness (mm)
Phibo Cr-Co copings	0,030	0,050	Yes
Phibo Cr-Co 2-7 copings	0,030	0,070	Yes
Phibo Cr-Co 8-14 copings	0,030	0,100	Yes
Phibo Ti copings	0,030	0,050	Yes
Phibo Ti 2-7 copings	0,030	0,070	Yes
Phibo Ti 8-14 copings	0,030	0,100	Yes
Phibo Zr 1-3 copings	0,020	0,090	Yes
Phibo Zr 4-7 copings	0,020	0,120	Yes
Phibo Wax copings	0,030	0,100	Yes
Phibo Wax 2-7 copings	0,030	0,070	Yes
Phibo Wax 8-14 copings	0,030	0,090	Yes
Phibo PMMA 2-7 copings	0,030	0,120	Yes
Phibo PMMA 8-14 copings	0,030	0,050	Yes
Phibo PMMA 2-7 copings	0,030	0,070	Yes
Phibo PMMA 8-14 copings	0,030	0,100	Yes
E.Max Cad HT Coping Ant	0,030	0,050	Yes
E.Max Cad LT Coping Post	0,030	0,050	Yes
E.Max Cad LT Coping Ant	0,030	0,050	Yes
E.Max Cad MO Coping Ant	0,030	0,050	Yes
E.Max Cad MO Coping Post	0,030	0,050	Yes

Coping coverage			
Nombre	Edge line adjustment (mm)	Angle adjustment (deg)	Wall thickness (mm)
Phibo Cr-Co copings	0,200	65,000	0,450
Phibo Cr-Co 2-7 copings	0,200	65,000	0,450
Phibo Cr-Co 8-14 copings	0,200	65,000	0,450
Phibo Ti copings	0,200	65,000	0,450
Phibo Ti 2-7 copings	0,200	65,000	0,450
Phibo Ti 8-14 copings	0,200	65,000	0,450

Coping coverage			
Name	Edge line adjustment (mm)	Angle adjustment (deg)	Wall thickness (mm)
Phibo Zr 1-3 copings	0,200	7,000	0,500
Phibo Zr 4-7 copings	0,200	70,000	0,500
Phibo Zr 8-14 copings	0,200	70,000	0,500
Phibo Wax copings	0,200	65,000	0,450
Phibo Wax 2-7 copings	0,200	65,000	0,450
Phibo Wax 8-14 copings	0,200	65,000	0,450
Phibo PMMA 2-7 copings	0,200	65,000	0,450
Phibo PMMA 8-14 copings	0,200	65,000	0,450
E.Max Cad HT Coping Ant	0,200	65,000	0,850
E.Max Cad LT Coping Post	0,200	65,000	1,000
E.Max Cad LT Coping Ant	0,200	65,000	0,850
E.Max Cad LT Coping Post	0,200	65,000	1,000
E.Max Cad MO Coping Ant	0,200	65,000	0,850
E.Max Cad MO Coping Post	0,200	65,000	1,000

CEMENTED SPACE	ADDITIONAL CEMENTED SPACE	EDGE LINE ADJUSTMENT	ANGLE ADJUSTMENT	LAYER THICKNESS

The number of pieces and/or intermediate structure connections in Phibo® CAD-CAM is unlimited for Titanium, Cobalt-Chrome, PMMA and wax for all types of structures. In the case of the Zirconia, it is unlimited for cement retained structures and limited to a maximum of seven pieces in screw retained structures. Only one cemented piece is possible with the IPS e-max® CAD system.

The customized abutments (for any material) must have a maximum height of 18 mm. In some special cases the maximum height may be 14 mm (these cases are indicated in Dental SystemTM modelling software; the software indicates the limitation when designing the abutment).

This package leaflet was revised and approved on 2016-09-19.

INSUSOCADCAM\_rev005

### ESPAÑOL

### ESTRUCTURAS INTERMEDIAS PHIBO® CAD-CAM INFORMACION IMPORTANTE. LEA DETENIDAMENTE