

phibo^φ

GROUP

AUREA EVO CATALOG

2025




VERSION: 00
CODE: PM-00001



AUREA EVO SYSTEM

phibo^φ
GROUP

RECOMMENDATIONS FOR INSERTION SECTORS





Mismatch:
0.25 mm

Ø 3.5 mm



Mismatch:
0.28 mm

Ø 4.3 mm



Mismatch:
0.53 mm

Ø 4.8mm



Mismatch:
0.88 mm

Ø 5.5 mm

NP




RP



RP



WP



AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Aurea Evo Narrow Platform Implant Ø 3.5mm Length	8.5mm	EVO NP 085
	10.0mm	EVO NP 100
	11.5mm	EVO NP 115
	13.0mm	EVO NP 130
	14.5mm	EVO NP 145



DESCRIPTION		REFERENCE
Aurea Evo Regular Platform Implant Ø 4.3mm Length	8.5mm	EVO RP 085
	10.0mm	EVO RP 100
	11.5mm	EVO RP 115
	13.0mm	EVO RP 130
	14.5mm	EVO RP 145



DESCRIPTION		REFERENCE
Aurea Evo Regular Platform Implant Ø 4.8mm Length	8.5mm	EVO RP 48085
	10.0mm	EVO RP 48100
	11.5mm	EVO RP 48115
	13.0mm	EVO RP 48130
	14.5mm	EVO RP 48145



DESCRIPTION		REFERENCE
Aurea Evo Wide Platform Implant Ø 5.5mm Length	8.5mm	EVO WP 085
	10.0mm	EVO WP 100
	11.5mm	EVO WP 115
	13.0mm	EVO WP 130



DESCRIPTION		REFERENCE
Aurea Evo Narrow Platform Healing Abutment	3mm	EVO NP 01.3
	4mm	EVO NP 01.4
	5mm	EVO NP 01.5
	6mm	EVO NP 01.6



DESCRIPTION		REFERENCE
Aurea Evo Regular Platform Healing Abutment	3mm	EVO RP 01.3
	4mm	EVO RP 01.4
	5mm	EVO RP 01.5
	6mm	EVO RP 01.6



DESCRIPTION		REFERENCE
Aurea Evo Wide Platform Healing Abutment	3mm	EVO WP 01.3
	4mm	EVO WP 01.4
	5mm	EVO WP 01.5
	6mm	EVO WP 01.6




DESCRIPTION		REFERENCE
Aurea Evo Narrow Platform Implant Freasable Abutment	1mm	EVO NP 20.1
	2mm	EVO NP 20.2
	3mm	EVO NP 20.3



DESCRIPTION		REFERENCE
Aurea Evo Regular Platform Implant Freasable Abutment	1mm	EVO RP 20.1
	2mm	EVO RP 20.2
	3mm	EVO RP 20.3



DESCRIPTION		REFERENCE
Aurea Evo Wide Platform Implant Freasable Abutment	1mm	EVO WP 20.1
	2mm	EVO WP 20.2
	3mm	EVO WP 20.3



DESCRIPTION		REFERENCE
Aurea Evo Narrow Platform 17° freasable Angled Abutment	1.5mm	EVO NP 23.0
	2.5mm	EVO NP 23.1



DESCRIPTION		REFERENCE
Aurea Evo Regular Platform 17° freasable Angled Abutment	1.5mm	EVO RP 23.0
	2.5mm	EVO RP 23.1
Aurea Evo Regular Platform 30° freasable Angled Abutment	1.5mm	EVO RP 63.0
	2.5mm	EVO RP 63.1

AUREA EVO REFERENCES

DESCRIPTION		REFERENCE
Aurea Evo Non-Hexed Temporary Abutment	NP	EVO NP 53.0
	RP	EVO RP 53.0
	WP	EVO WP 53.0
Aurea Evo Temporary Hexed Abutment	NP	EVO NP 53.1
	RP	EVO RP 53.1
	WP	EVO WP 53.1



DESCRIPTION		REFERENCE
Aurea Evo Implant Analog	NP	EVO NP 03.0
	RP	EVO RP 03.0
	RP 4.8	EVO RP48 03.0
	WP	EVO WP 03.0

AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Implant Impression Coping Aurea Evo O.T	NP	EVO NP 02.0
	RP	EVO RP 02.0
	WP	EVO WP 02.0
Implant Impression Coping Aurea Evo C.T.	NP	EVO NP 02.2
	RP - WP	EVO RW 02.2

DESCRIPTION		REFERENCE
Sinergia Library AUREA EVO (DI)	NP	166S
	RP	167S
	RP 48	219S
	WP	168S



AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Titanium Base Axis® PHIBO AUREA EVO Non-Hexed	NP	INTAEVONPR
	RP - WP	INTAEVORWR
Titanium Base Axis® PHIBO AUREA EVO Hexed	NP	INTAEVONP
	RP - WP	INTAEVORW

DESCRIPTION		
Titanium Base Axis® PHIBO AUREA EVO (P) Non-Hexed	NP	INTAEVOPNPR
	RP	INTAEVOPRPR
	WP	INTAEVOPWPR
Titanium Base Axis® PHIBO AUREA EVO (P) Hexed	NP	INTAEVOPNP
	RP	INTAEVOPRP
	WP	INTAEVOPWP



DESCRIPTION		REFERENCE
Titanium Base Axis® PHIBO AUREA EVO ANGULADO (P) Non-Hexed	RP - WP	INTAEVOANR
Titanium Base Axis® PHIBO AUREA EVO ANGULADO (P) Hexed	RP - WP	INTAEVOAN


DESCRIPTION		REFERENCE
Screw for Phibo Axis® abutment connection	M1,4	TORPPIM14
	M1,6	TOREXPIM16
	M1,8	TOREXPIM18



AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Aurea Evo Narrow Platform Abutment	1mm	EVO NP 11.1
	2mm	EVO NP 11.2
	3mm	EVO NP 11.3
	4mm	EVO NP 11.4
Aurea Evo Regular Platform Abutment	1mm	EVO RP 11.1
	2mm	EVO RP 11.2
	3mm	EVO RP 11.3
	4mm	EVO RP 11.4
Aurea Evo Wide Platform Abutment	1mm	EVO WP 11.1
	2mm	EVO WP 11.2
	3mm	EVO WP 11.3
	4mm	EVO WP 11.4




DESCRIPTION		REFERENCE
Aurea Evo Narrow Platform 17° Angled Abutment	1.5mm	EVO NP 22.0
	2.5mm	EVO NP 22.1
Aurea Evo Regular Platform 17° Angled Abutment	1.5mm	EVO RP 22.0
	2.5mm	EVO RP 22.1
Aurea Evo Regular Platform 30° Angled Abutment	1.5mm	EVO RP 62.0
	2.5mm	EVO RP 62.1

AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Aurea Evo Short Healing Cap	NP	EVO NP 49.1
	RP	EVO RP 49.1
	WP	EVO WP 49.1
Aurea Evo Short Angled Healing Cap	NP-RP	EVO NR 30.1



DESCRIPTION		REFERENCE
Hexed Aurea Evo Abutment Direct Coping O.T.	NP	EVO NP 26.1
	RP	EVO RP 26.1
	WP	EVO WP 26.1
Aurea Evo Hexed Angled Abutment Metallic Direct Coping O.T.	NP-RP	EVO NR 27.1



DESCRIPTION		REFERENCE
Sinergia Library AUREA EVO	NP	169S
	RP	170S
	WP	171S
Sinergia Library AUREA EVO (Angled Abutment)	NP-RP-WP	172S

AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Aurea Evo Non-Hexed Temporary Titanium Cap	NP	EVO NP 21.1
	RP	EVO RP 21.1
	WP	EVO WP 21.1
Aurea Evo Hexed Temporary Titanium Cap	NP	EVO NP 12.1
	RP	EVO RP 12.1
	WP	EVO WP 12.1



DESCRIPTION		REFERENCE
Aurea Evo Hexed Temporary Angled Ti Abutment Cap	NP - RP	EVO NR 24.1
Aurea Evo No-Hexed Temporary Angled Ti Abutment Cap	NP - RP	EVO NR 24.0



DESCRIPTION		REFERENCE
Aurea Evo Abutment Analog	NP	EVO NP 18.0
	RP	EVO RP 18.0
	WP	EVO WP 18.0
Aurea Evo Angled Abutment Analog	NP-RP	EVO NR 28.0



AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Aurea Evo Hexed Casting Cylinder	NP	EVO NP 14.1
	RP	EVO RP 14.1
	WP	EVO WP 14.1
Aurea Evo Non-Hexed Casting Cylinder	NP	EVO NP 14.2
	RP	EVO RP 14.2
	WP	EVO WP 14.2



DESCRIPTION		REFERENCE
Aurea Evo Angled Abutment Hexed Casting Cylinder	NP - RP	EVO NR 48.1
Aurea Evo Angled Abutment Non-Hexed Casting Cylinder	NP - RP	EVO NR 48.2



DESCRIPTION		REFERENCE
Aurea Evo Platform Waxing Screw	NP-RP-WP	EVO NW 19.0



DESCRIPTION		REFERENCE
Aurea Evo Abutment Prosthetic Screw	NP-RP-WP	EVO NW 15.0



DESCRIPTION		REFERENCE
CLICK & FIX Abutment for AUREA EVO NP	2mm	2017 C
	3mm	2018 C
	4mm	2019 C
	5mm	2020 C
CLICK & FIX Abutment for AUREA EVO RP	2mm	2023 C
	3mm	2024 C
	4mm	2025 C
	5mm	2026 C



DESCRIPTION		REFERENCE
CLICK & FIX Retention kit	1 unit	8519 C
CLICK & FIX Titanium Hole	10 units	CAZTI
CLICK & FIX Black Retention	10 units	RETNG
CLICK & FIX Blue Retention	10 units	RETAZ
CLICK & FIX Pink Retention	10 units	RETRS
CLICK & FIX Red Retention	10 units	RETRJ
CLICK & FIX Transparent Retention	10 units	RETTR
CLICK & FIX Separator	10 units	SEPSI

AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Precision Drill	STD	EVO 23000



DESCRIPTION		REFERENCE
Pilot Drill	Ø2.0mm	EVO 20000



DESCRIPTION		REFERENCE
Surgical Drill Aurea Evo	NP Ø3.1mm	EVO 31000
	NP Ø3.4mm	EVO 34000
	RP Ø3.8mm	EVO 38000
	RP Ø4.1mm	EVO 41000
	RP Ø4.4mm	EVO 44000
	RP Ø4.6mm	EVO 46000
	WP Ø5.2mm	EVO 52052
	WP Ø5.4mm	EVO 54000

AUREA EVO REFERENCES




DESCRIPTION		REFERENCE
Depth Indicator Drill	Ø2.0 mm	EVO 00200



DESCRIPTION		REFERENCE
Torque Wrench	STD	172.0172



DESCRIPTION		REFERENCE
Open end Wrench	STD	172.0001



DESCRIPTION		REFERENCE
Aurea Evo Surgical Box	STD	EVO 00001



DESCRIPTION		REFERENCE
Drill stop Aurea Evo	NP 8.5mm	TOP NP 085
	NP 10.0mm	TOP NP 100
	NP 11.5mm	TOP NP 115
	NP 13.0mm	TOP NP 130
	NP 14.5mm	TOP NP 145
	RP 8.5mm	TOP RP 085
	RP 10.0mm	TOP RP 100
	RP 11.5mm	TOP RP 115
	RP 13.0mm	TOP RP 130
	RP 14.5mm	TOP RP 145
	RP48 8.5mm	TOP RP48 085
	RP48 10.0mm	TOP RP48 100
	RP48 11.5mm	TOP RP48 115
	RP48 13.0mm	TOP RP48 130
	RP48 14.5mm	TOP RP48 145
	WP 8.5mm	TOP WP 085
	WP 10.0mm	TOP WP 100
	WP 11.5mm	TOP WP 115
	WP 13.0mm	TOP WP 130

AUREA EVO REFERENCES




DESCRIPTION		REFERENCE
Bone Tap Aurea Evo	NP	EVO 01035
	RP	EVO 01043
	RP 4.8	EVO 01048
	WP	EVO 01055



DESCRIPTION		REFERENCE
Circular Scalpel	Ø3.70	152.0001
	Ø4.70	152.0002
	Ø6.00	152.0003



DESCRIPTION		REFERENCE
Ratchet Adapter,- to Mounter	Short	172.0100
	Long	172.0300



DESCRIPTION		REFERENCE
Manual Hex Driver Abutment ConHex TSH®	Short	174.2301
	Medium	174.2302



DESCRIPTION		REFERENCE
Contra-Angle Adapter	Short	173.0100
	Long	173.0300

AUREA EVO REFERENCES



DESCRIPTION		REFERENCE
Contra-Angle Hex Driver 1.25mm	Short	173.1251
	Medium	173.1252



DESCRIPTION		REFERENCE
Manual Hex Driver 1.25mm	Short	174.1251
	Medium	174.1252
	Long	174.1253



DESCRIPTION		REFERENCE
Axis® Hexalobular Angled Tool - to Wrench	Short	172.1103
	Medium	172.1104
	Long	172.1105
Contra-Angle Axis® Hexalobular Angled Tool	Medium	173.1106



DESCRIPTION		REFERENCE
Hexalobular tool - to wrench	Short	172.1101
	Medium	172.1102



SYSTEM AUREA EVO

phibo^φ
GROUP